

BINOCULAR ACCOMMODATION  
AUTO REF/KERATO METER  
**WAM-5500**  
ADVANCED

### Specification

Refraction Measurement	Sphere -22~+22D (VD=0)
	-30~+17.4D (VD=12) (0.01/0.12/0.25DStep)
	Cylinder 0~±10D (0.01/0.12/0.25DStep)
	Axis angle 0~180° (1°Step)
	Vertex Distance 0, 10, 12, 13.5, 15mm
Measurement of Corneal Radius	Minimum Pupil Diameter $\phi$ 2.3mm
	Corneal radius 5.0~10.0mm (0.01mmStep)
	Refractive power 33.75~67.5D (0.01/0.12/0.25DStep)
	Cylindrical power 0~±9D
Pupil Diameter Measurement	Axis angle 0~180° (1°Step)
Pupillary Distance	$\phi$ 2~ $\phi$ 8mm (0.1mmStep)
Printer	Measurement Range 85mm (1mmStep)
Internal Monitor	Thermal printer with Automatic Cutter (Width 57mm)
Movable Distance	5.6 inch LCD Display (Color)
Movable Distance of Chinrest	Back/Force $\pm$ 17mm Right/Left $\pm$ 43mm Up/Down $\pm$ 15mm
Overall Dimension	$\pm$ 30mm
Weight	(W) 327mm $\times$ (D) 496mm $\times$ (H) 515mm
Output	About 20kg
Rated Voltage	RS232C Interface
Consumption	Video Output (NTSC)
Power Save	100~240V
	50/60Hz
	80VA
	OFF, 3, 5, 10 min. (Selectable)

- Specifications and design are subject to change without prior notice for improvement.
- The part of screen is composed photograph.
- The color of instrument on the catalogue and the real product might be different.
- We judge that the LCD Monitor is qualified if the total "lit pixel" and/or "missing pixel" is less than five (excluding non-sharp one and less than half one).
- Medical equipment manufacturer license No.37B3X00002

**Rexxam** CE 0197  
Rexxam Co., Ltd.  
958 Ikeuchi, Konan-cho, Takamatsu, Kagawa 761-1494 Japan



EC REP

**REXXAM ELECTRONICS IRELAND CO.LIMITED**  
Donore Road, Drogheda, County Louth, Ireland  
Tel: +353-41-9839700 Fax: +353-41-9839702

**SHIGIYA**  
**SHIGIYA MACHINERY WORKS LTD.**

GS Division <http://www.GrandSeiko.com>

5378 Minoshima, Fukuyama, Hiroshima, 721-8575 Japan  
Tel: +81-84-981-5735 Fax: +81-84-953-6758


### Accessories

- Test eye / 1 pc ● Power cable / 1 pc ● Printer paper / 3 rolls
- Fuse 2A / 2 pcs ● Chinrest paper / 1 pack (1000 sheets)
- Chinrest pin / 2 pcs ● Dust cover / 1 pc ● Fixation target / 1 pc
- Fogging lens / 1 pc ● Occluder / 1 pc
- Near point measurement unit / 1 set ● Contact lens holder / 1 pc
- WCS-1 HI SPEED MODE Control software

### Option

- WMT-1 Moving target system ● WSM-1 Status monitor
- MDC-1 The measurement data collection software
- Magnet stopper ● Electric table

### Overall Dimension



Hardware requirement  
OS : Japanese/English Version of Microsoft Windows Vista/7 /8  
CPU : Core2 or better (recommended)  
Memory : 1GB or more (recommended)  
(Mentioned company and product names are the trademark or the registered trademark of each company.)



For the proper and safety use, please read the Instruction Manual thoroughly before using the product.

● Agency

Binocular, Dynamic, Refractive, Pupil Diameter Measurements  
More improvements are now available

**NEW**

BINOCULAR ACCOMMODATION  
AUTO REF/KERATO METER

**WAM-5500**

ADVANCED



  
**Grand Seiko.com**



# ADVANCED ACCOMMODATION MEASUREMENT

## Binocular, Dynamic, Refractive, Pupil Diameter Measurements

### New CPU with the faster High Speed Mode measurement From 5 times /second to 6 times /second

BINOCULAR ACCOMMODATION  
AUTO REF/KERATO METER  
**WAM-5500**  
ADVANCED

In addition to the Binocular Refractor and Keratometer measurement functions, the WAM-5500 can measure the Refraction and the Pupil Diameter simultaneously. These are advanced features provided on the earlier WR-5100K AutoRef/Keratometer version. With this function, it can be used to confirm the change on the accommodation function before and after the surgery. In addition, it can be used to confirm presbyopia and the accommodation error due to fatigue as a result of VDT work or other eye fatigue instances. By connecting with the computer, you can measure refractive value at a rapid, continuous speed, thus allowing you to measure the accommodation response dynamically.



Pupil Diameter  
Selection of A.M.MODE  
Pupil Image

No. 00001	
NAME 2014 1 1 12:00	
VD=0	
<R> FAR	Selection of A.M.MODE
-4.75 -0.25 62	Ref Value
-4.75 0.00 5.0	Pupil Diameter
-4.75 0.00 5.0	
-4.75 0.00 5.0	
-4.75 0.00 5.0	Representative Ref Value
-4.75 0.00 5.0	Average Pupil Diameter
<L> -2.0D/50cm	Selection of A.M.MODE
-4.50 -0.75 90	Ref Value
-4.50 -0.75 90	Pupil Diameter
-4.50 -0.75 90	
-4.50 -0.75 90	
-2.5D/40cm	Selection of A.M.MODE
-4.50 -0.75 89	(When MODE is changed)
-4.50 -0.75 90	
-3.0 D/33cm	Selection of A.M.MODE
-4.50 -0.75 90	(When MODE is changed)
-4.50 -0.75 90	
-4.50 -0.75 90	
PD = 64	
Grand Seiko.com	
WAM-5500 S/N:000265	

Left Eye Measurement Data Right Eye Measurement Data

#### A.M. MODE

It measures the pupil diameter at 0.1mm step (minimum 2mm) simultaneously with refractive value. The measurements are then presented on the monitor. Subsequently, you can print the results and check the constriction of pupil. Studies have indicated that pupil constriction is evident when the eye is accommodating. Thus, the pupil diameter measurement function allows you to more precisely ascertain whether the subject's eyes are accommodating or not. Another option is to present the image via video thus giving you precise image analysis of the pupil constriction as well. The ideal refractive value for 50, 40, 33, 25, and 20cm is displayed, thus you can compare the patient data with the ideal value more easily as it relates to the actual distance of the Near Point Card.




#### Easier measurement of IOL patient

Compare to before, unnecessary reflection from IOL is reduced and it leads to more assured measurement of IOL Patient.

#### Easier access to the eye ground image

Able to confirm the eye ground image right after the measurement (Freeze Image function). It is useful to check the reason of non-measurement and confirm the good or bad measurement.


In the case IOL is not at the proper position, it sometimes cannot be measured. With this function, the alignment can be adjusted easily for the patient that cannot be measured even with the IOL mode.



#### HI SPEED MODE

By connecting with the computer and using our dedicated software - WCS-1, you can measure the SE value and Pupil Diameter at 0.16 seconds per step. In addition, you can analyze this data with spreadsheet software like Microsoft Excel.

Unlike the previous model, WR-5100K, this Hi Speed Mode allows you to more dynamically measure the eye. Additionally, in this Mode, since the refractive value is calculated on the entire meridian instead of a specific one, you can understand the refractive value precisely.



● Example ; 33 years old male  
This is the data of moving the Near Point Card every 20 seconds. With this data, the subject eyes are accommodating from the ideal value of -2.00D to -5.00D. Additionally, the pupil diameter is constricting from 4.0mm to 3.2mm. (Blinking at the same interval.)

#### Error handling mechanism during the Hi Speed Mode is improved

Up till now, the measurement could not be taken for a few seconds after the patient blinks and the difference would be too big compare to the previous ones. With the improved algorithm, it would judge the blink and move to the next measurement immediately without losing the data for a few seconds.

#### Measurement with Both Eyes

With both eyes open, the natural condition, the measurement is taken. Thus, there would be less Instrumental Myopia, and allows you to achieve more precise measurement.

#### Easy Measurement with Open View Window

The Open View Window, which adopted the Half Mirror, makes it possible to take a measurement without the subject being aware as there is no distraction of view. For the infant/children, it gives less pressure, and you can measure easily and quickly. Additionally, it allows you to go through the measurement process smoothly since you can observe the condition of a subject directly.

#### Any Target

You can select from any target, as you like. Needless to say, it can be an eye-test chart. Or for the infant/children, you can choose from what they are interested in, such as toy and cartoon. It helps the measurement process go more smoothly. Additionally, since you can put a target at any distance, you can measure middle and near target as well.

#### Indication of SE (Spherical Equivalent)

You can confirm the change of accommodation as the SE is indicated on the monitor each time (3 data).

#### Easy Over-Refracton Measurement

Since it is adopting the Open View Window, you can measure a subject with glasses and contact lens as well as a subject with the IOL, easily. It also allows you to confirm the prescription.



Right eye data	Date & time
No. 00001	
NAME 2014 1 1 12:00	
VD=12	
<R> SPH CYL AX	Refractive data
-4.75 0.00	
-4.75 0.00	
-4.75 0.00	
-4.75 0.00	Representative value
<R> mm D AX	Corneal data
R1 7.62 44.25	
R2 7.62 44.25	
AVE 7.62 44.25	
CYL 0.00	
<L> SPH CYL AX	Refractive data
-4.75 0.00	
-4.75 0.00	
-4.75 0.00	
-4.75 0.00	Representative value
<L> mm D AX	Corneal data
L1 7.62 44.25	
L2 7.62 44.25	
AVE 7.62 44.25	
CYL 0.00	
PD = 67 NPD=62	Pupillary distance
Grand Seiko.com	
WAM-5500 S/N:000265	
Left eye data	