

CXL UNBOUND™

THE C-EYE™ DEVICE

BY EMAGINE



Slit Lamp Mount

For modern CXL: fits on Haag Streit and Zeiss



Pulsed and Continuous Light

For various keratoconus & infectious keratitis settings.



Table Mount

For classic CXL: table mount accessory available.



Superior Battery Technology

Built to last for thousands of treatments.



Range of Intensities

Delivers 3, 9, 18, 30 and 36 mW/cm².



Thickness-Adjusted


Beam profile delivers more energy in the corneal periphery.




Customized settings

Or choose from a number of preset clinical protocols.

ONE DEVICE TO PERFORM BOTH...



CXL at the slit lamp
The modern way



CXL in the operating room
The classic way

Innovative

EMAGine's C-Eye™ device brings cross-linking technology to the slit lamp. Simple, efficient and safe from the novice to the expert user.

Versatile

Perform CXL at the slit lamp. Or use the C-Eye™ device just like any other CXL device: on a table mount in your operating room.

Portable

Work in different places? Take the C-Eye™ device wherever you go. Charge it via USB-C with your smartphone or laptop charger.



RIBO-KER™

ONE RIBOFLAVIN FOR KERATOCONUS AND KERATITIS TREATMENT



EMAGine is proud to offer Ribo-Ker™, its exclusive CE-marked riboflavin.

No more confusion with different types of riboflavin: EMAGINE has developed one solution for all applications

- ✓ 0.1% hypoosmolar riboflavin
- ✓ CE marked
- ✓ Ready for epi-off, epi-on and iontophoresis applications
- ✓ Can be used on normal thickness and thin corneas

EMAGine AG
Gubelstrasse 12
6300 Zug, Switzerland
T: +41 41 560 91 92 / F: +41 41 560 91 99

The CE Mark for the C-Eye™ device is pending and is not currently available in the US.

 **EMAGINE**
www.emagine-eye.com